INSTRUCTIONS FOR USE

SUPERFLEX® EAGLE FLEXPLUG® EAGLEPLUG® TAPERED-SHAFT™ FLOW CONTROLLER SUPEREAGLE® PLUG 1™ EAGLEPLUG® TEARFLOW™

STERILE R   2797

DESCRIPTION

The SuperFlex®, Eagle FlexPlug®, EaglePlug®, Tapered-Shaft™ Flow Controller, SuperEagle®, PLUG 1™ and EaglePlug® TearFlow™ are designed for insertion into the punctal opening to block or modulate tear drainage through the canaliculus. They are molded of a soft medical grade silicone, with a colorant of titanium dioxide, and come pre-loaded on a disposable, inserter/dilator for ease of insertion.

The patented designs feature a flattened rim which remains outside the punctum and a central shaft which tapers outward toward a pointed nose. EagleVision® Tapered-Shaft™ Flow Controllers, and EaglePlug® TearFlow™ have an open nose. The Eagle FlexPlugs have a ribbed shaft. The tapered shaft is designed to create a vector force to keep the plug snugly in the punctal opening.

These products are intended for use by a physician practicing in ophthalmology.

The following instructions apply to the SuperFlex®, Eagle FlexPlug®, EaglePlug®, Tapered-Shaft™ Flow Controller, SuperEagle®, Plug 1™ and EaglePlug® TearFlow™.

CUSTOM FITTING

The EagleVision Punctal Gauging Systems are designed to size the punctal opening, ensuring a custom fit with maximum retention and patient comfort. Gauge sizes are the exact diameter stated on the gauge. Plug sizes are the same diameters as the gauge, with a slightly oversized nose for maximum plug retention.

GAUGING (SIZING)

Note: *Determine the proper plug size prior to dilation.*

Place patient in a reclining position for best results. Anesthetize the area around the punctum using a topical anesthetic. Begin with the smallest size gauge. Gently insert, taking care not to go past the neck of the gauge (*fig.1*).

The gauge head should slightly flex the punctal ring. If the smallest gauge enters the punctal ring with substantial resistance, use a smaller plug. If there is no resistance, use a larger size gauge until it moderately flexes the punctal ring when entering

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
(fig. 1)

and exiting the punctal opening. Use a plug the same size as the gauge. If the largest gauge enters the punctal ring with little or no resistance, use a larger plug. If there is sizing uncertainty, always select

a plug one size smaller. Remember, the plug nose is slightly larger than the gauge. (If patient is to receive plugs in both eyes, gauge both eyes due to possible variances in size.)

DILATION

(fig. 2)

Gently insert the dilator end of the inserter/dilator instrument no more than 1-2mm into the vertical canaliculus (*fig.2*). Remove and immediately insert the punctum plug per the instructions

below. If the plug is not inserted within 15 to 30 seconds, the dilation may have to be repeated.

INSERTION

It is recommended to lubricate the plug on the pre-loaded inserter pin with saline solution or artificial tears to ease insertion. Hold the instrument as you would a pencil, between the thumb and middle finger, with your index finger just above the release button. DO NOT DEPRESS THE BUTTON AT THIS TIME.

Ensure the plug is fully seated on the inserter wire by pressing lightly with a blunt object on the nose of the plug. It may be necessary to adjust the trigger so that the length of exposed wire between the plug and the inserter tip is no more than 0.5mm, aiding in proper placement of the punctum plug.

Using a gentle, downward rotational motion, insert the plug until the rim is flush against the punctal opening (*fig. 3, 4*). Once the plug is properly positioned, fully depress and hold the release button on the inserter. The inserter pin has now been withdrawn

from the plug. Remove the disposable instrument while holding the button down (*fig.5*). Dispose of used inserter in “sharps” container.

PLUG REMOVAL

Gently grasp the plug, under the rim, as far down the shaft as possible with suitable, small forceps. Gently pull the plug from the punctal opening. Do not try to remove by grasping the edges of the exposed rim.

INDICATIONS FOR USE

The SuperFlex®, Eagle FlexPlug®, EaglePlug®, SuperEagle® and Plug 1™ are designed for punctal and canalicular occlusion in patients with chronic dry eye. Tapered-Shaft™ Flow Controller and EaglePlug® TearFlow™ are designed to regulate the flow of lacrimal fluid through the punctum and canaliculus.

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Punctum plugs may be used in the treatment of dry eye syndrome and the dry eye components of various ocular surface diseases such as corneal ulcers, conjunctivitis, pterygium, blepharitis, keratitis, red lid margins, recurrent chalazions, recurrent corneal erosion, filamentary keratitis.

In addition, Punctum Plugs, Tapered-Shaft™ Flow Controller, or EaglePlug® TearFlow™ may be useful in treating various dry eye related problems by increasing comfort for LASIK eye surgery or contact lens induced dry eye and retention/enhancement of ocular medications on the eye.

CONTRAINDICATIONS

The SuperFlex®, Eagle FlexPlug®, EaglePlug®, Tapered-Shaft™ Flow Controller, SuperEagle®, PLUG 1™ and EaglePlug® TearFlow™ are contraindicated for patients with known silicone sensitivity, infective conjunctivitis, dacryocystitis, inflammation of the eyelid, irritated or infected eyes or epiphora.

CAUTIONS AND POSSIBLE COMPLICATIONS

- U.S Federal law restricts this device for sale by or on the order of a licensed medical practitioner.
- If injecting an anesthetic agent in the region of the canaliculus, maintain a minimum 5mm distance between the injection path and the angular vessels.
- Do not dilate the punctal opening more than 2mm to avoid rupture of the punctum or canaliculus. Rupture may cause pain, increase the risk of infection or cause the plug to fall out due to an enlarged punctal opening. If rupture occurs, delay the insertion of the plug until the wound heals.
- Consider plug removal if the patient experiences irritation, epiphora, erythema, pruritus, swelling, infection or pain that persists longer than several days. Conjunctiva abrasion is more commonly produced by upper lid plugs.
- If protrusion or partial extrusion occurs (more common in patients over 60 years of age), reset the plug. If condition persists, consider replacing with a larger size plug. Total extrusion with loss of the plug may occur if the punctum is over dilated or from the patient wiping the eye.
- Intrusion of the plug into the canaliculus can occur from excessive insertion force or, in rare cases, occur spontaneously. Suppurative canaliculitis or pyogenic granuloma may result.
- Use of silicone punctum plugs may enhance the effect of ocular medication used on the eye. Dosage may need to be reduced accordingly, depending on the type of medication used.
- REGULATORY INFORMATION - Any serious incident that occurs in relation to this product should be reported to the manufacturer, and to the healthcare authority of the country in which the incident occurred.

WARNINGS

- Supplied Sterile.
- Sterility not guaranteed if package has been opened or damaged.

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- DO NOT RESTERILIZE OR REUSE. Eagle Vision assumes no liability for products which have been resterilized by health care facilities. The device may be compromised if it is reused due to damage when explanted or degradation due to resterilization. There is also a potential of infection.
- Dispose of used inserter in “sharps” container.
- For patient safety, do not eliminate spacing between the plug and tip of the inserter; the stainless steel wire could penetrate the plug and result in injury to the punctum.
- Do not force punctal plug into punctum if plug appears too large to fit. The pin may move forward, perforate the highly compressed plug, and harm patients punctal opening.

STORAGE CONDITIONS

Store at Room Temperature.

Can withstand -20°C to 55°C.

EC REP

Advena Ltd, Tower Business Centre, 2nd Flr, Tower Street, Swatara, BKR 4013, Malta

	- Do Not Re-use
	- “Use by” (expiration) date
	- Do not use if package is damaged
STERILE R	- Sterilized using irradiation
	- Single sterile barrier system
LOT	- Batch code
REF	- Catalog Number
	- Consult instructions for use
MD	- Medical device
	- Manufacturer
	- Date of manufacture
Rx Only	- For use by or on the order of a physician
EC REP	- Authorized European Representative

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©2023 Katena Products, Inc. 930030K Rev G 2023-04-30