



blink SINGLE-USE
INSTRUMENTS

Single-Use Ophthalmic Instruments

OPHTHALMOLOGY

 **corzamedical**

For ophthalmic surgeons and healthcare facilities seeking the highest standards in patient safety and surgical precision, Corza Medical provides an extensive range of Blink branded ophthalmic single-use instruments.





Blink Medical has become widely recognised as a market-leading company in the manufacture and supply of single-use ophthalmic instruments and custom packs.

Blink single-use instruments are sold globally by Corza Medical, alongside premium reusable ophthalmic instrument brands such as Katena and ASICO.

We provide our customers with a multitude of solutions to meet their ever-changing needs.

Table of Contents

Product Group	Page Number
Tying Forceps	4
General Forceps	5
Fixation Forceps	6
Jewellers Forceps • Mosquito Forceps	7
Capsulorhexis Forceps	8
Sprung Scissors	9
Speculums	10 - 13
Iris and Tenotomy Scissors	13
Manipulators • Choppers	14
Double-Ended Manipulators/Choppers	15
Needle Holders	16
Refractive • Corneal	17
Markers	18
Sub Tenon Anaesthesia	19
Oculoplastics	20 - 21
Standard Packs	22
General	22


Tying Forceps

Straight

Tying Forceps

10-9008 6mm Tying Platform,
O/L 110mm

10-9009 5mm Tying Platform,
O/L 90mm



Angled

Tying Forceps

11-6015 5mm Platform Angled,
O/L 105mm


11-6016 5mm Platform Angled,
O/L 85mm



Curved

Tying Forceps

11-6012 5mm Platform Angled,
O/L 108mm



General Forceps

Adson Forceps

Fine Toothed Tip, O/L 120mm

10-1305



Adson Brown Forceps

O/L 120mm

10-1307



Max Fine Forceps Curved

Smooth Pointed Tip, O/L 108mm

11-6009



Moorfields Forceps

16mm Serrated Platform, O/L 102mm

10-1301



Moorfields Suture Utility Forceps

14mm Narrow Serrated Platform, O/L 115mm

10-1311



Epilation Forceps

Angled Tip, O/L 91mm

10-1300




Fixation Forceps

Hoskins Forceps Straight

10-8008 0.3mm Notched, O/L 90mm

10-8009 0.3mm Notched, O/L 108mm



Pierse Colibri Forceps

0.3mm Notched, O/L 84mm

11-7009



Bonn Forceps Straight

10-1106 0.25mm Toothed, O/L 106mm

10-1107 0.25mm Toothed, O/L 90mm

10-1108 0.12mm Toothed, O/L 102mm



Colibri Forceps

(108mm version shown)

11-7010 0.25mm Toothed, O/L 84mm

11-7011 0.25mm Toothed, O/L 108mm



Lister Utility Forceps

0.25mm Toothed, O/L 112mm

10-1302



St Martins Forceps

10-1203 0.6mm Toothed, O/L 90mm

10-1204 0.3mm Toothed, O/L 90mm



Jewellers Forceps

Jewellers Forceps #1

O/L 120mm

10-2001



Jewellers Forceps #3

O/L 120mm

10-2003



Jewellers Forceps #5

O/L 112mm

10-2005



Jewellers Forceps #5 Angled

O/L 108mm

11-2005



Jewellers Forceps #7

O/L 117mm

11-2007



Mosquito Forceps

Mosquito Forceps

10-1006 Straight,
O/L 100mm

11-1005 Curved,
O/L 100mm



10-1006




11-1005



Capsulorhexis Forceps


Capsulorhexis Forceps Utrata

- II-5012** Round Handle, Straight Shaft, Cystotome Tip, O/L 105mm
- II-5013** Round Handle, Curved Shaft, Cystotome Tip, O/L 105mm
- II-5014** Round Handle, Curved Shaft, Grasping Tip, O/L 105mm
- II-5015** Round Handle, Straight Shaft, Grasping Tip, O/L 105mm



Capsulorhexis Forceps Utrata

- II-5050** Flat Handle, Straight Shaft, Cystotome Tip, O/L 105mm
- II-5051** Flat Handle, Curved Shaft, Cystotome Tip, O/L 105mm
- II-5052** Flat Handle, Straight Shaft, Grasping Tip, O/L 105mm
- II-5053** Flat Handle, Curved Shaft, Grasping Tip, O/L 105mm



Capsulorhexis Forceps Inamura

- II-5021** Flat Handle, Curved Shaft, Cystotome Tip, O/L 102mm
- II-5011** Flat Handle, Straight Shaft, Cystotome Tip, O/L 102mm



Capsulorhexis Forceps Inamura

- II-5020** Round Handle, Curved Shaft, Cystotome Tip, O/L 105mm



Sprung Scissors

Westcott Scissors

- 20-4013** Straight, O/L 115mm
- 20-4015** Straight, Sharp Tip, O/L 115mm
- 21-4010** Curved, O/L 115mm
- 21-4014** Curved, Sharp Tip, O/L 115mm
- 21-4011** Curved, Dolphin Nose, O/L 115mm
- 21-5001** Curved, O/L 100mm



Curved Tip



Straight Tip



Dolphin Nose



Vannas Scissors

- 20-6000** Straight, O/L 90mm
- 21-6001** Curved, O/L 90mm



20-6000



21-6001



Ong Scissors

- 21-3010** Angled Tip 9mm, O/L 108mm




Speculums

Neonatal Speculum

Angled Temporal,
Width 20mm,
B/L 7.4mm,
O/L 35mm


31-3100



Neonatal Speculum with finger-loops

Width 12mm,
B/L 5.5mm and 3.5 mm,
O/L 32mm


31-3103



Paediatric Speculum

Width 23mm,
B/L 7mm,
O/L 40mm

31-3102



Paediatric Speculum Spring

Solid Blade,
Width 26mm,
B/L 7.4mm,
O/L 35mm

30-3002



Speculums

Clarks Speculum

Solid Blade, Temporal,
Width 10-60mm,
B/L 13.5mm,
O/L 95mm

31-3004



Burratto Speculum

Width 15-30mm,
B/L 14mm,
O/L 40mm

30-3001



Cook Speculum

Solid Blade,
Width 13-40mm,
B/L 13mm,
O/L 40mm

31-3003



Solid Plastic Bladed Speculum

Width 38mm,
B/L 14.5mm,
O/L 40mm

31-4001




Speculums

Barraquer Closed Loop

30-2004 Straight, Width 30mm,
B/L 15mm,
Length 43mm

31-2003 Temporal, Width 30mm,
B/L 15mm,
Length 42mm

31-2002 Nasal, Width 30mm,
B/L 15mm,
Length 42mm



Kratz-Barraquer Open Loop

31-1006 Temporal, Width 30mm,
B/L 13.5mm,
Length 40mm

30-1007 Straight, Width 30mm,
B/L 13.5mm,
Length 42mm

31-1005 Nasal, Width 30mm,
B/L 13.5mm,
Length 40mm

31-1009 Temporal, Curved, Width 30mm
B/L 13.5mm
Length 42mm



31-1009



31-1006

Leibinger Speculum

Wire Loop, Temporal,
Width 15-50mm,
B/L 14mm,
O/L 85mm

31-3002



Speculums

Leibinger Speculum

Solid Blade, Temporal,
Width 15-50mm,
B/L 16mm,
O/L 80mm

31-3000



Barraquer Speculum

Solid Blade, Temporal,
Width 30mm,
B/L 13mm,
O/L 44mm

30-4000



Iris and Tenotomy Scissors

Iris Scissors

20-7001 Iris, Straight, Blunt Tip,
O/L 100mm

20-7002 Iris, Straight, Sharp Tip,
O/L 100mm


21-7004 Iris, Curved, Sharp Tip,
O/L 100mm



Tenotomy Scissors

21-7005 Tenotomy, Curved,
Pointed Tip, O/L 115mm


21-7006 Tenotomy, Curved,
Blunt Tip, O/L 115mm



Manipulators


Manipulators

- 40-5002** Sinskey Lens Hook, Straight
- 41-5001** Sinskey Lens Hook, Angled
- 41-2005** Nucleus Y Shaped, Angled Horizontal Axis
- 41-2006** Nucleus Y Shaped, Angled Vertical Axis
- 41-3004** Spatula Manipulator
- 41-4014** Kuglen Hook, Angled
- 41-8009** Mushroom Tipped Manipulator
- 41-1007** IOL Lester Manipulator
- 41-3005** Spatula Neuhann Manipulator
- 41-4013** Connor Wand



Choppers

- 41-6007** Drysdale Paddle Chopper
- 41-4005** Sadiq Chopper
- 41-4008** Koch Chopper, 1.4mm
- 41-4010** Rosen Chopper, 1.2mm
- 41-4012** Seibel Chopper
- 41-4009** Nagahara Chopper
- 41-7008** Greenes Chopper
- 41-3006** McPhearson Chopper



Double-Ended Manipulators and Choppers

Sadiq Chopper/Mushroom Tipped Manipulator
41-9012



Sinskey Manipulator/Mushroom Tipped Manipulator
41-9010



Rosen Chopper/Mushroom Tipped Manipulator
41-9013



Nagahara Chopper/Connor Wand
41-9014



Needle Holders

Needle Holder Castroviejo

12mm Tying Platform, Straight With Lock,
O/L 135mm

90-7001



Needle Holder Castroviejo

12mm Tying Platform, Straight,
O/L 135mm

90-7002



Needle Holder Castroviejo

12mm Tying Platform, Delicate Curved,
O/L 135mm

91-7003



Needle Holder Troutman

6mm Tying Platform, Curved,
O/L 105mm

91-7004



Refractive and Corneal

Refractive Toric Alignment Marker

3, 6, 9 O'Clock

91-6001



Refractive Toric Degree Gauge

Bevelled face with markings in 5° increments

11.7mm Internal diameter

91-6000



Refractive Toric IOL Marker

11mm Wide

91-6002



Epithelial Spatula

7mm flat blade length

51-0001



Epithelial Spatula

Rounded Tip

51-0003



Alcohol Well

9.5mm Diameter

51-0002



Stutman Spud and Burr

O/L 126mm

90-1000



Epley Scleral Depressor

O/L 136mm

51-0004



Markers

Double-ended Marker

Flat Handle, 3.5mm-4mm, O/L 113mm

90-1015



Braunstein Marker

Round Handle, 3.5mm-4mm, O/L 80mm

90-1017



Caliper & Rule

0-20mm Measurement Span

90-1014



Intravitreoretinal Caliper/Fixator

Round Handle, 3.5mm, O/L 70mm

90-1016



Sub Tenon Anaesthesia

Conjunctival Probe

O/L 83mm


CP01



CP Sub Tenon Set

Barraquer Speculum, Angled Temporal - 31-2006,
Moorfields Forceps - 10-1301,
Conjunctival Probe - CP01
Supplied in box of 10

270



Sub Tenon Set Fine Tip

Westcott Scissors - 21-4010,
Tip Speculum Kratz-Barraquer, Angled Temporal - 31-1006,
Moorfields Forceps - 10-1301
Supplied in boxes of 8

002



Oculoplastics

Liebreich Probes

- 90-2007** 0000+000 O/L 125mm
- 90-2008** 000+00 | O/L 125mm
- 90-2009** 00+0 | O/L 125mm
- 90-2010** 1+2 | O/L 125mm
- 90-2011** 3+4 | O/L 125mm



Curette Double-ended

2mm & 3mm Cups, O/L 130mm

90-3007



Curette Double-ended

2mm & 5mm Cups, O/L 130mm


90-3008



Curette Double-ended

2mm & 4mm Cups, O/L 130mm

90-3009



Clamp Chalazion

- 90-4005** Small, 11mm Diameter Clamp, O/L 90mm
- 90-4006** Medium, 18mm x 11mm Clamp, O/L 80mm
- 90-4007** Large, 23mm x 15mm Clamp, O/L 90mm



Desmares Entropian Clamp

Clamp Diameter,
24mm x 16mm, O/L 95mm

- 90-4008** Right
- 90-4009** Left



Oculoplastics

Meibomian Gland Expression Forceps

3mm Tip Width, O/L 105mm

11-3015



Desmarres Retractor

14mm Wide Blade, O/L 130mm

90-9013



Graefe Strabismus Hook

10mm Tip Length, O/L 140mm

90-5004



Alexandra Muscle Hook

11mm Tip Length, O/L 130mm

90-5005



Dilator Punctal Double-Ended

Tapered Tips, O/L 138mm

90-1011



Dilator Nettle Ship

20mm Tapered Tip, O/L 118mm

90-1012



Bulldog Clamp

O/L 52mm

90-8012



Standard Packs

Chalazion Set Small-Clamp

Chalazion Clamp Small - 90-4005,
Curette Double Ended Small - 90-3007,
No.11 Scalpel,

092



Chalazion Set Medium-Clamp

Chalazion Clamp Medium - 90-4006,
Curette Double Ended Medium - 90-3009,
No.11 Scalpel,

III



General

Thornton Ring

13mm Diameter

90-1020



Backhaus Towel Clamp

Opens to 8cm, O/L 95mm

90-5050





Order Online: corzaeye.com

Order by Phone: 800-225-1195

For more information, visit corzaeye.com

Distributed by Corza Ophthalmology, (Katena Products Inc)

Corza Ophthalmology (Katena Products, Inc.), 6 Campus Drive, Suite 310, Parsippany, NJ 07054, USA
Phone: 973.989.1600, Toll free (USA only): 800.225.1195, Fax: 973.989.8175

All trademarks are owned by the company and its affiliates.

©2025 Corza Medical All Rights Reserved. COR-CATBLINK-062124-Rev1