

● ● ● DIAMOND
● ● ● BI-ASPHERIC LENSES



katena®



DIAMOND BI-ASPHERIC LENSES

Katena's Diamond series of reusable lenses feature a patented optical design at an extraordinary value. Each lens possesses high performance multi layer anti-reflective and hard coatings. An ultralight metallic alloy with silicone grip makes each lens ergonomically friendly. These lenses offer exceptional small pupil performance and low image distortion.



K31-2350 90D LENS

This reusable 90D lens provides an image magnification of 0.64x and a field of view of 75°.



K31-2365 28D LENS

This reusable 28D lens provides an image magnification of 2.1x and a field of view of 53°.



K31-2355 78D LENS

This reusable 78D lens provides an image magnification of 0.77x and a field of view of 80°.



K31-2370 20D LENS

This reusable 20D lens provides an image magnification of 3.0x and a field of view of 53°.



K31-2360 60D LENS

This reusable 60D lens provides an image magnification of 0.96x and a field of view of 68°.