



**Sharpoint** MICROSURGICAL  
KNIVES

**Single-Use Ophthalmic  
Surgical Knives**

 **corzamedical**

# Setting a new standard in knives.

## Single-Use Ophthalmic Surgical Knives

Since 1969, Sharpoint has focused on one objective: redefining what surgeons expect of their microsurgical knives. Built on decades of knife engineering expertise, Sharpoint knives are designed not only for exceptional sharpness, but for the complete performance experience—how the knife engages tissue, how efficiently it initiates an incision, how smoothly it advances, and how reliably it performs from case to case.

Your operative experience, control and consistency in every case is the result of four key factors—sharpness, penetration force, geometry, and drag, which we call Blade Acuity™. At Corza Medical, we are committed to applying our legacy knowledge and experience through the surgeon-driven design and our industry-leading manufacturing processes that deliver a knife as consistent as the hand that holds it.



# The Factors of Blade Acuity™

Your operative experience is influenced by more than the sharpness of your blade. Control and consistency in every case requires an optimized balance of key factors that together we call Blade Acuity™.

## 1 Sharpness

Fineness of the cutting edge itself, typically measured at the microscopic level.  
The knife has to be sharp every time.

## 2 Penetration Force

Initial force needed to begin a cut, influenced by tip design and edge configuration.  
A predictable entry—no surprise resistance, no over-correction.

## 3 Geometry

Blade shape, edge angle, and overall profile affect how it interacts with anatomy.  
Edge angle and profile held to spec, not to chance.

## 4 Drag

Resistance encountered as the blade moves through material, affected by surface finish, blade thickness, and edge geometry.  
The same smooth glide through tissue, blade after blade.

# Table of Contents

|                                            |     |
|--------------------------------------------|-----|
| Safety Knives. . . . .                     | 3   |
| Slit/Keratome Knives . . . . .             | 4-5 |
| ClearTrap™ Trapezoid Knives . . . . .      | 6   |
| Implant Knives . . . . .                   | 6   |
| Paracentesis/MVR/Sideport Knives . . . . . | 7   |
| Vitrectomy Knives . . . . .                | 8   |
| Crescent Knives . . . . .                  | 8   |
| Groove & LRI Knives. . . . .               | 9   |
| Specialty Knives & Blades . . . . .        | 9   |

# Safety Knives

## Safety Slit

SINGLE  
BEVEL

5 per box



| Item #   | Size    | Angle | Depth Indicator |
|----------|---------|-------|-----------------|
| 62-2260G | 2.2 mm  | 35°   | Yes             |
| 62-2460G | 2.4 mm  | 35°   | Yes             |
| 62-2560G | 2.5 mm  | 35°   | Yes             |
| 62-2760G | 2.75 mm | 35°   | Yes             |
| 62-2860G | 2.8 mm  | 35°   | Yes             |
| 62-3060G | 3.0 mm  | 35°   | Yes             |

## Safety Slit—CC3D+Advanced

DOUBLE  
BEVEL

5 per box




| Item #   | Size   | Angle | Depth Indicator |
|----------|--------|-------|-----------------|
| 62-2250G | 2.2 mm | 35°   | Yes             |
| 62-2450G | 2.4 mm | 35°   | Yes             |

## Safety Sideport

MVR

5 per box



| Item #   | Size   | Angle    |
|----------|--------|----------|
| 68-2010G | 1.0 mm | Straight |

## Safety Slit

DOUBLE  
BEVEL

5 per box




| Item #   | Size    | Angle | Depth Indicator |
|----------|---------|-------|-----------------|
| 62-2442G | 2.4 mm  | 35°   | Yes             |
| 62-2542G | 2.5 mm  | 35°   | Yes             |
| 62-2742G | 2.75 mm | 35°   | Yes             |
| 62-3042G | 3.0 mm  | 35°   | Yes             |

## Safety Sideport Stab

STAB

5 per box



| Item #   | Angle |
|----------|-------|
| 62-1501G | 15°   |

## Safety Crescent

BEVEL-UP

5 per box



| Item #   | Size   | Angle        |
|----------|--------|--------------|
| 64-1030G | 2.3 mm | 55° Bevel-Up |

# Slit/Keratome Knives

## Slit/Keratome

BEVEL-UP

6 per box



| Item #  | Size    | Angle | Depth Indicator | Blade/Tip |
|---------|---------|-------|-----------------|-----------|
| 72-2261 | 2.2 mm  | 40°   | Yes             | —         |
| 72-2461 | 2.4 mm  | 40°   | Yes             | —         |
| 72-2561 | 2.5 mm  | 40°   | Yes             | —         |
| 72-2661 | 2.65 mm | 40°   | Yes             | —         |
| 72-2761 | 2.75 mm | 40°   | Yes             | 72-2762   |
| 72-2861 | 2.85 mm | 40°   | Yes             | 72-2862   |
| 72-2961 | 3.0 mm  | 40°   | Yes             | —         |
| 72-3261 | 3.2 mm  | 40°   | Yes             | —         |

## Clear Corneal Knives

DOUBLE BEVEL

6 per box



| Item #  | Size    | Angle | Depth Indicator |
|---------|---------|-------|-----------------|
| 72-1361 | 1.3 mm  | 40°   | Yes             |
| 72-2263 | 2.2 mm  | 40°   | Yes             |
| 72-2463 | 2.4 mm  | 40°   | Yes             |
| 72-2663 | 2.65 mm | 40°   | Yes             |
| 72-2763 | 2.75 mm | 40°   | Yes             |
| 72-2863 | 2.85 mm | 40°   | Yes             |
| 72-3061 | 3.0 mm  | 40°   | Yes             |
| 72-4963 | 4.9 mm  | 40°   | Yes             |

## Slit Knives

SINGLE BEVEL

6 per box



| Item #  | Size    | Straight |
|---------|---------|----------|
| 72-1001 | 1.0 mm  | Straight |
| 72-2801 | 2.85 mm | Straight |
| 72-2901 | 3.0 mm  | Straight |

## Slit/Keratome

DOUBLE BEVEL

6 per box



| Item #  | Size    | Angle | Depth Indicator |
|---------|---------|-------|-----------------|
| 72-1831 | 1.8 mm  | 40°   | Yes             |
| 72-2231 | 2.2 mm  | 40°   | Yes             |
| 72-2431 | 2.4 mm  | 40°   | Yes             |
| 72-2531 | 2.5 mm  | 40°   | Yes             |
| 72-2631 | 2.65 mm | 40°   | Yes             |
| 72-2831 | 2.8 mm  | 40°   | Yes             |
| 72-2931 | 3.0 mm  | 40°   | Yes             |
| 72-3231 | 3.2 mm  | 40°   | Yes             |

# Slit/Keratome Knives *(continued)*

## CC3D+ Slit/Keratome

CC3D+ 6 per box



| Item #  | Size    | Angle | Depth Indicator |
|---------|---------|-------|-----------------|
| 72-1650 | 1.6 mm  | 40°   | Yes             |
| 72-1850 | 1.8 mm  | 40°   | Yes             |
| 72-2050 | 2.0 mm  | 40°   | Yes             |
| 72-2250 | 2.2 mm  | 40°   | Yes             |
| 72-2450 | 2.4 mm  | 40°   | Yes             |
| 72-2550 | 2.5 mm  | 40°   | Yes             |
| 72-2650 | 2.65 mm | 40°   | Yes             |
| 72-2750 | 2.75 mm | 40°   | Yes             |
| 72-2850 | 2.85 mm | 40°   | Yes             |
| 72-3050 | 3.0 mm  | 40°   | Yes             |
| 72-3250 | 3.2 mm  | 40°   | Yes             |

## Clear Corneal—Parallel Sides

DOUBLE BEVEL 6 per box



| Item #  | Size    | Angle | Depth Indicator |
|---------|---------|-------|-----------------|
| 72-2635 | 2.65 mm | 40°   | Yes             |
| 72-2735 | 2.75 mm | 40°   | Yes             |
| 72-2835 | 2.85 mm | 40°   | Yes             |
| 72-3035 | 3.0 mm  | 40°   | Yes             |

## Slit Knives—Parallel Sides

BEVEL-UP 6 per box



| Item #  | Size   | Angle | Depth Indicator | Blade/Tip |
|---------|--------|-------|-----------------|-----------|
| 72-2912 | 3.0 mm | 40°   | No              | 72-2762   |
| 72-3212 | 3.2 mm | 40°   | No              | 72-2862   |

# ClearTrap™ Trapezoid Knives

## ClearTrap™

TRAPEZOID 6 per box



| Item #  | Size         | Angle |
|---------|--------------|-------|
| 75-1214 | 1.2–1.4 mm   | 40°   |
| 75-1617 | 1.5–1.7 mm   | 40°   |
| 75-1620 | 1.5–2.0 mm   | 40°   |
| 75-1720 | 1.7–2.0 mm   | 40°   |
| 75-2022 | 2.0–2.2 mm   | 40°   |
| 75-2522 | 2.5–2.9 mm   | 40°   |
| 75-2723 | 2.75–3.0 mm  | 40°   |
| 75-2724 | 2.75–3.2 mm  | 40°   |
| 75-2725 | 2.75–3.5 mm  | 40°   |
| 75-2824 | 2.85–3.2 mm  | 40°   |
| 75-2826 | 2.85–3.75 mm | 40°   |

## IQ Geometry ClearTrap™

BEVEL-UP 6 per box



| Item #  | Size           | Angle |
|---------|----------------|-------|
| 77-1725 | 1.7–2.2–2.5 mm | 40°   |

# Implant Knives

## Pilot Tip Implant Knives

BEVEL-UP 6 per box



| Item #  | Size   | Angle |
|---------|--------|-------|
| 76-5361 | 3.2 mm | 40°   |
| 76-5161 | 5.2 mm | 40°   |
| 76-5561 | 5.5 mm | 40°   |

## Implant Knives

BEVEL-UP 6 per box



| Item #  | Size   | Angle |
|---------|--------|-------|
| 76-4061 | 4.0 mm | 40°   |
| 76-5261 | 5.2 mm | 40°   |

# Paracentesis/MVR/Sideport Knives

## ClearPort™ Paracentesis 6 per box



| Item #  | Size   | Angle | Straight/<br>Angle |
|---------|--------|-------|--------------------|
| 78-2310 | 0.8 mm | 45°   | Angled             |
| 78-2020 | 1.0 mm | 45°   | Angled             |
| 78-2010 | 1.1 mm | 45°   | Angled             |
| 78-2019 | 1.4 mm | 45°   | Angled             |

## Stab Knives 6 per box

RESTRICTED  
DEPTH



| Item #  | Angle | Depth  | Blade/Tip |
|---------|-------|--------|-----------|
| 72-1531 | 15°   | 3.0 mm | 72-1532   |
| 72-1551 | 15°   | 5.0 mm | 72-1552   |
| 72-3036 | 30°   | 3.5 mm | —         |

## Stab Knives 6 per box



| Item #  | Size  | Config.      | Blade/Tip |
|---------|-------|--------------|-----------|
| 72-1501 | 15°   | Straight     | 72-1502   |
| 72-1591 | 15°   | Right-Angled | —         |
| 72-2201 | 22.5° | Straight     | 72-2202   |
| 72-3001 | 30°   | Straight     | 72-3002   |
| 72-4501 | 45°   | Straight     | 72-4502   |

# Vitrectomy Knives

## Vitrectomy

6 per box



| Item #  | Gauge    | Straight/Angle |
|---------|----------|----------------|
| 71-1901 | 19 Gauge | Straight       |
| 71-2001 | 20 Gauge | Straight       |
| 71-2301 | 23 Gauge | Straight       |
| 71-2501 | 25 Gauge | Straight       |
| 71-2061 | 20 Gauge | Angled 40°     |

# Crescent Knives

## Sharptome Crescent Knives

6 per box



| Item #  | Size   | Straight/<br>Angled | Bevel        | Blade/<br>Tip |
|---------|--------|---------------------|--------------|---------------|
| 74-1000 | 2.0 mm | Straight            | Bevel-Up     | 74-5000       |
| 74-1010 | 2.0 mm | Angled              | Bevel-Up     | 74-5062       |
| 74-1031 | 2.0 mm | Angled              | Double-Bevel | —             |

## Mini-Crescent Knives

6 per box



| Item #  | Size    | Straight/<br>Angled | Bevel        |
|---------|---------|---------------------|--------------|
| 74-6910 | 1.25 mm | Angled              | Double-Bevel |

# Groove & LRI Knives

## LRI Knives

PRECISION  
DEPTH

6 per box



| Item #  | Depth  |
|---------|--------|
| 72-0303 | 300 µm |
| 72-3503 | 350 µm |
| 72-5503 | 550 µm |
| 72-6003 | 600 µm |

## Restricted Depth Knives

RESTRICTED

6 per box



| Item #  | Depth  |
|---------|--------|
| 72-0301 | 300 µm |
| 72-0401 | 450 µm |

# Specialty Knives & Blades

## Spoon Knives

DOUBLE  
BEVEL

6 per box



| Item #  | Size   | Straight/<br>Angled | Bevel        |
|---------|--------|---------------------|--------------|
| 74-2261 | 2.2 mm | Angled              | Double-Bevel |
| 74-3061 | 3.0 mm | Angled              | Double-Bevel |

## Blade & Tip Handle

1 per box



| Item #  | Length | Type      |
|---------|--------|-----------|
| 92-9003 | 16 mm  | Universal |

## Specialty Blades

6 per box



| Item #  | Size           | Straight/<br>Angled | Description                    |
|---------|----------------|---------------------|--------------------------------|
| 78-6810 | Needle Blade   | Straight            | 4.0 mm cutting edge            |
| 78-5700 | Blade          | Straight            | Specialty Blade/Tip            |
| 78-6600 | Lamellar Blade | Angled              | Bevel-Up, 2.3 mm               |
| 78-6900 | Micro-Blade    | Straight            | Double-Bevel, Sharp all around |



78-6810



78-5700



78-6600



78-6900

## Notes

This image shows a single sheet of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.



## USA

Order Online: [corzaeye.com](https://corzaeye.com)

Order by Phone: 800-225-1195

## UK

[salesuk@corza.com](mailto:salesuk@corza.com)

Order by Phone: 0121 386 8433

## International

For all international inquiries:

[inquiries.corzaeye@corza.com](mailto:inquiries.corzaeye@corza.com)

**For more information, visit [corzaeye.com](https://corzaeye.com)**

Distributed by Corza Ophthalmology, (Katena Products Inc)

Corza Ophthalmology (Katena Products, Inc.), 6 Campus Drive, Suite 310, Parsippany, NJ 07054, USA  
Phone: 973.989.1600, Toll free (USA only): 800.225.1195, Fax: 973.989.8175

Manufactured by Surgical Specialties México S. DE R.L. DE C.V.  
Corredor Tijuana-Rosarito 2000, #24702-B, Ejido Francisco Villa, Tijuana, B.C., C.P. 22235  
All trademarks are owned by the company and its affiliates.

©2026 Corza Medical. All Rights Reserved. COR-SHARPOINT-CAT-0726 R0