

# Tutoplast®



## The Patch Graft Standard

### Tutoplast® Processed Allograft Advantages:

- Appropriate for use with all glaucoma drainage devices
- Allows for room temperature storage
- Five year shelf life
- Biocompatible

### Tutoplast® Pericardium

- Nominal thickness of 400 microns
- Multi-directional matrix for superior surgical handling and suture utility
- Low profile prevents graft bulkiness over the implant

### Tutoplast® Sclera

- High profile graft delays conjunctival erosion
- Multi-directional matrix for superior surgical handling and suture utility
- Pre-cut sizes for various applications

### Tutoplast® Clear Cornea

- Native clear cornea with adjacent, peripheral sclera to form a sturdy, customizable, cosmetically appealing option

## Applications (Refer to product labeling for complete instructions for use)

### Glaucoma

- Drainage implant/valve implantation
- Repair of leaking and/or overfiltering bleb
- Repair corneo-scleral fistulas

### Retina

- Repair exposed scleral buckles

### Cataract

- Phaco “burns”

### Oculoplastics

- Eyelid reconstruction
- Cover of exposed orbital implants
- Gold eyelid weight cover

**More than 5 million implants have been sterilized through the Tutoplast® process with 0 confirmed incidence of implant-associated infection.**

# Ordering Information

## Tutoplast® Pericardium

| Description   | Item # |
|---------------|--------|
| 1.0 x 1.0 cm  | 68251  |
| 1.5 x 1.5 cm  | 68250  |
| 2.0 x 3.0 cm  | 68260  |
| 1.0 x 2.5 cm  | 68256  |
| 4.0 x 5.0 cm  | 68252  |
| 5.0 x 5.0 cm  | 68257  |
| 2.0 x 5.0 cm  | 68255  |
| 6.0 x 12.0 cm | 68254  |

## Tutoplast® Sclera

| Description  | Item # |
|--------------|--------|
| 0.5 x 0.8 cm | 68333  |
| 0.6 x 1.0 cm | 68337  |
| 1.0 x 3.0 cm | 68334  |
| Whole Sclera | 68332  |

## Tutoplast® Clear Cornea

| Description            | Item # |
|------------------------|--------|
| 1/4 Cornea with Sclera | 68328  |
| 1/3 Cornea with Sclera | 68329  |
| 1/2 Cornea with Sclera | 68330  |

Reimbursement:

Coding for reimbursement is dependent upon the indication and geographic region.

Please use the Katena reimbursement hotline for further guidance: 888.700.9005

For more information, or to place an order, contact your local representative.

Tutoplast® is a registered trademark of Tutogen Medical GmbH



800-225-1195 • [www.katena.com](http://www.katena.com)

