

RETINA 200

SINGLE USE LENS FOR RETINAL TREATMENT



K30-1330

FEATURES

The Retina 200 represents a breakthrough in lens optics. The optical design of the lens permits imaging and treatment of the far peripheral retina with less distortion.

- Expansive field of view to observe and treat the entire retina
- Observation of the retina including the ora serrata
- Single use lens provided in a sterile pouch
- Anti-reflection coating
- Supplied in quantities of 10

APPLICATIONS

- Facilitate the diagnosis and treatment of diabetic eye disease
- Assist in the diagnosis of age related macular degeneration
- Pan retinal photocoagulation

SPECIFICATIONS

Image magnification:	0.60x
Laser spot magnification:	1.67x
Field of view	200°