


SMALL 3 MIRROR LENS

SINGLE USE LENS FOR DIAGNOSIS & LASER



K30-1002



SIZE COMPARISON

FEATURES

The small 3 mirror lens offers eye care providers a choice of single use 3 mirror lenses without sacrificing the exceptional optics of Katena's comprehensive lens selection.

- Excellent all-purpose lens especially sized for smaller palpebral fissure & pediatric patients
- Mirror optimized for viewing through small pupils
- Excellent choice for patients difficult to dilate
- Single use lens provided in a sterile pouch
- Supplied in quantities of 10

APPLICATIONS

- General purpose lens with a central optic to view posterior pole, and mirrors to view & treat the anterior chamber angle, peripheral, and arcades regions

SPECIFICATIONS

Image magnification:	1.00x
Laser spot magnification:	1.00x
Central fundus view	30°