

The Vidaurri Fluid Retention Ring

for **PRK & CXL**



K20-2135 Vidaurri Fluid Retention Ring,
8.7mm diameter

Hands free solution control!


- Retain solution on cornea alone
- Minimize leakage with double walled suction ring
- Reduce the amount of solution used in every case
- Eliminate cross linking effects on limbal and conjuntival cells
- Supplied sterile and disposable

Developed with Jesus Vidaurri, MD of Monterrey, Mexico

 **katena**
DESIGNED FOR SIGHT™



cross-hairs for centration



solution retained on cornea