


# Hands down, a better gonio lens than you have ever handled.

The NEW Hands-Free Gonio from Corza Ophthalmology.




**Convenient Lens Tab**  
Allows ease of lens positioning to desired view.

**Glass-like Optics Material**  
Sharp, clear, detailed vision of the surgical area.

**Lightweight and Compact**  
Designed to minimize pressure on the globe.

**No Silicone Flange**  
Designed for greater instrument access.



The new Hands-Free Gonio provides exceptional clarity with direct views of angle structures, enhanced hands-free design characteristics and unique lightweight features.

The Hands-Free Gonio delivers precise lens handling during intraoperative bimanual gonioscopy procedures and is the clear choice for optimizing procedure efficiencies — all as a convenient, single-use disposable lens.

Contact your local sales representative or visit [www.corzaeye.com](http://www.corzaeye.com).

# Optimizes intraoperative control and visibility.

## Hands-Free Gonio Lens



### LENS SPECIFICATIONS

<b>Gonio Mag</b>	1.2x	• Small, lightweight and hands-free, designed to deliver outstanding balance and secure positioning with minimal pressure on globe
<b>Lens Weight</b>	0.5g	
<b>Lens Height</b>	8mm	• Innovative glass-like optic offers a 1.2x magnified view of angle structures and an overall wide view of the anterior chamber
<b>Lens Width</b>	11mm	
<b>Lens Length</b>	11mm	• No flange lens designed to provide greater instrument access and bimanual flexibility
<b>Item#</b>	K30-1407	• Hands-free lens optimizes surgical efficiency putting focus on the procedure and other surgical tasks
<b>Pack</b>	10 per box	• Single-use, sterile device, conveniently packaged and ready for use, with no need for cleaning, sterilization, or disinfection